

Product Datasheet

BAG2 antibody (orb19039)

Description	Goat polyclonal antibody to BAG2
Species/Host	Goat
Reactivity	Bovine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	BAG2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 25kDa band observed in Human HeLa lysates (calculated MW of 23.7kDa according to NP_004273.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	CSKTLQQAESRFN
Clonality	Polyclonal
MW	23.7
Entrez	9532
Dilution Range	ELISA: 1:8000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Pig Testis lysa...



Western blot analysis of HeLa lysate usi...