

Product Datasheet

Human GABAA coupled receptor Antibody (orb1905787)

Description	Human GABAA coupled receptor Antibody
Species/Host	Camelus
Tested Applications	IF, IHC
Target	Human GABAA coupled receptor
Preservatives	25 mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, 0.09 % proclin
Form/Appearance	Liquid/ lyophilized
Concentration	1mg/mL
Storage	Storage Condition: 4°C; Do not freeze
Note	For research use only
Isotype	Camelid VHH
Clonality	Monoclonal
Purity	> 95% as determined by SDS-PAGE
MW	15 kDa
Expiration Date	12 months from date of receipt.