

Product Datasheet

Human transmembrane protein 30A Antibody (orb1905793)

Description	Human transmembrane protein 30A Antibody
Species/Host	Camelus
Tested Applications	IF, IHC
Target	Human transmembrane protein 30A
Preservatives	25 mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, 0.09 % proclin
Form/Appearance	Liquid/ lyophilized
Concentration	1mg/mL
Storage	Storage Condition: 4°C; Do not freeze
Note	For research use only
Isotype	Camelid VHH
Clonality	Monoclonal
Purity	> 95% as determined by SDS- PAGE
MW	15 kDa
Expiration Date	12 months from date of receipt.