

## Product Datasheet

# Human transmembrane protein 30A Antibody (orb1905796)

<b>Description</b>	Human transmembrane protein 30A Antibody
<b>Species/Host</b>	Camelus
<b>Tested Applications</b>	IF, IHC
<b>Target</b>	Human transmembrane protein 30A
<b>Preservatives</b>	25 mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, 0.09 % proclin
<b>Form/Appearance</b>	Liquid/ lyophilized
<b>Concentration</b>	1mg/mL
<b>Storage</b>	Storage Condition: 4°C; Do not freeze
<b>Note</b>	For research use only
<b>Isotype</b>	Camelid VHH
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95% as determined by SDS- PAGE
<b>MW</b>	15 kDa
<b>Expiration Date</b>	12 months from date of receipt.