

Product Datasheet

human Tau protein Antibody (orb1905802)

| | |
|----------------------------|---|
| Description | human Tau protein Antibody |
| Species/Host | Camelus |
| Conjugation | Unconjugated |
| Tested Applications | IF, IHC |
| Target | human Tau protein |
| Preservatives | 25 mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, 0.09 % proclin |
| Form/Appearance | Liquid/ lyophilized |
| Concentration | 1mg/mL |
| Storage | Storage Condition: 4°C; Do not freeze |
| Note | For research use only |
| Isotype | Camelid VHH |
| Clonality | Monoclonal |
| Purity | > 95% as determined by SDS-PAGE |
| MW | 15 kDa |
| Expiration Date | 12 months from date of receipt. |