

Product Datasheet

anti-EGFR VHH antibody (orb1905877)

Description	anti-EGFR VHH antibody
Species/Host	Camelus
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Target	Human Epidermal Growth Factor Receptor 1 (HER1/EGFR)
Preservatives	25 mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, 0.09 % proclin
Form/Appearance	Liquid
Concentration	1mg/mL
Storage	Storage Condition: 4°C; Do not freeze. Shelf Life: 1 year
Note	For research use only
Isotype	Camelid VHH
Clonality	Monoclonal
Purity	> 95% as determined by SDS-PAGE
MW	15 kDa
Expiration Date	12 months from date of receipt.