

Product Datasheet

Recombinant Rat Fractalkine/CX3CL1 (orb1905996)

Description	Recombinant Rat Fractalkine/CX3CL1
Endotoxins	Less than 1 EU/ug of rRtFractalkine/CX3CL1 as determined by LAL method.
Form/Appearance	Lyophilized from a 0.2 µm filtered concentrated solution in 2x PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.
Note	For research use only
Protein Sequence	QHLGMTKCN I TCHKMTSPIP VTL LIHYQLN QESCGKRAII LETRQHRHFC ADPKEKVVQD AMKHLDHQTA ALTRNG
Purity	> 95 % by SDS-PAGE and HPLC analyses.
Source	Escherichia coli
MW	Approximately 8.8 kDa, a single non-glycosylated polypeptide chain containing 76 amino acids and comprises only the chemokine domain of Rat Fractalkine.
Uniprot ID	O55145
Expiration Date	6 months from date of receipt.

