

Product Datasheet

Recombinant Rat VEGF Co-regulated Chemokine 1/CXCL17 (orb1905998)

Description	Recombinant Rat VEGF Co-regulated Chemokine 1/CXCL17
Endotoxins	Less than 1 EU/ug of rRtVCC-1/CXCL17 as determined by LAL method.
Form/Appearance	Lyophilized from a 0.2 µm filtered concentrated solution in 2 x PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.
Note	For research use only
Protein Sequence	SPNQEVARRH GDQHQAPRRW LWEGGQECDC KDWSLRVSKR KTTAVLEPPR KQPCDHVKG SEKKNRRQKH HRKSQRPSRT CQQFLKRCQL ASFTLPL
Purity	> 97 % by SDS-PAGE and HPLC analyses.
Source	Escherichia coli
MW	Approximately 11.5 kDa, a single non-glycosylated polypeptide chain containing 97 amino acids.
Expiration Date	6 months from date of receipt.