

## Product Datasheet

# Recombinant Rat gamma-Interferon Inducible Protein 10/CXCL10 (orb1906001)

|                         |  |
|-------------------------|--|
| <b>Description</b>      | Recombinant Rat gamma-Interferon Inducible Protein 10/CXCL10   |
| <b>Endotoxins</b>       | Less than 1 EU/ugof rRtIP-10/CXCL10 as determined by LAL method.   |
| <b>Form/Appearance</b>  | Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.   |
| <b>Storage</b>          | Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution. |
| <b>Note</b>             | For research use only  |
| <b>Protein Sequence</b> | IPLARTVRCT CIDFHEQPLR PRAIGKLEII<br>PASLSCPHVE IIATMKKNNE KRCLNPESEA<br>IKSLLKAVSQ RRSKRAP   |
| <b>Purity</b>           | > 95 % by SDS-PAGE and HPLC analyses.  |
| <b>Source</b>           | Escherichia coli   |
| <b>MW</b>               | Approximately 8.7 kDa, a single non-glycosylated polypeptide chain containing 77 amino acids.  |
| <b>Uniprot ID</b>       | <b>P48973</b>  |
| <b>Expiration Date</b>  | 6 months from date of receipt.   |

