

Product Datasheet

Recombinant Human Macrophage Inflammatory Protein-5/CCL15 (orb1906060)

Description	Recombinant Human Macrophage Inflammatory Protein-5/CCL15
Endotoxins	Less than 1 EU/ug of rHuMIP-5/CCL15 as determined by LAL method.
Form/Appearance	Lyophilized from a 0.2 um filtered concentrated solution in PBS.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.
Note	For research use only
Protein Sequence	QFINDAETEL MMSKLPLENP VVLNSFHFAA DCCTSYISQS IPCSLMKSYF ETSSECSKPG VIFLTKKGRQ VCAKPSGPGV QDCMKKLPY SI
Purity	> 97 % by SDS-PAGE and HPLC analyses.
Source	Escherichia coli
MW	Approximately 10.2 kDa, a single non-glycosylated polypeptide chain containing 92 amino acids.
Uniprot ID	Q16663
Expiration Date	6 months from date of receipt.

