

Product Datasheet

Recombinant Human Monocyte Chemotactic Protein-2/CCL8 (orb1906065)

Description	Recombinant Human Monocyte Chemotactic Protein-2/CCL8
Endotoxins	Less than 1 EU/ug of rHuMCP-2/CCL8 as determined by LAL method.
Form/Appearance	Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM PB, pH 7.4, 100 mM NaCl.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.
Note	For research use only
Protein Sequence	QPDSVSIPIT CCFNVINRKI PIQRLESYTR ITNIQCPKEA VIFKTKRGKE VCADPKERWV RDSMKHLDQI FQNLKP
Purity	> 96 % by SDS-PAGE and HPLC analyses.
Source	Escherichia coli
MW	Approximately 8.9 kDa, a single non-glycosylated polypeptide chain containing 76 amino acids.
Uniprot ID	P80075
Expiration Date	6 months from date of receipt.

