

Product Datasheet

Recombinant Human Interferon-alpha2b, Yeast (orb1906253)

Description	Recombinant Human Interferon-alpha2b, Yeast
Endotoxins	Less than 0.1 EU/ug of rHuIFN- α 2b, Yeast as determined by LAL method.
Form/Appearance	Lyophilized from a 0.2 μ m filtered solution in PBS, pH 7.4, with 0.02 % Tween-20.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 $^{\circ}$ C as supplied.- 1 month, 2 to 8 $^{\circ}$ C under sterile conditions after reconstitution.- 3 months, -20 to -70 $^{\circ}$ C under sterile conditions after reconstitution.
Note	For research use only
Protein Sequence	CDLPQTHSLG SRRTLMLLAQ MRRISLFSCL KDRHDFGFPQ EEFGNQFQKA ETIPVLHEMI QQIFNLFSTK DSSAAWDETL LDKFYTELYQ QLNDLEACVI QGVGVTTETPL MKEDSILAVR KYFQRITLYL KEKKYSPCAW EVVRAEIMRS FSLSTNLQES LRSKE
Purity	> 98 % by SDS-PAGE and HPLC analyses.
Source	Yeast
MW	Approximately 19.3 kDa, a single polypeptide chain containing 165 amino acids.
Uniprot ID	P01563
Expiration Date	6 months from date of receipt.

