

Product Datasheet

Biotinylated Human GAS6 (C-Avi-6His) (orb1906459)

Description	Biotinylated Recombinant Human Growth arrest-specific protein 6 is produced by our Mammalian express...
Form/Appearance	Supplied as a 0.2 µm filtered solution of 20mM PB, 8% Sucrose, 0.05% Tween 80, pH 7.0
Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-Avi-6His
Note	For research use only
Uniprot ID	Q14393
Expiration Date	6 months from date of receipt.