

Product Datasheet

Biotinylated Human Latent TGF-beta 3 (N-8His-Avi) (orb1906471)

Description	Biotinylated Recombinant Human Transforming Growth Factor Beta-3 Proprotein is produced by our Mamma...
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-8His-Avi
Note	For research use only
Uniprot ID	P10600
Expiration Date	6 months from date of receipt.