

Product Datasheet

Recombinant Human/Mouse/Rat Irisin/FNDC5 (C-6His) (orb1906492)

Description	Recombinant Human/Mouse/Rat Fibronectin type III domain-containing protein 5 is produced by our Mamm...
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Uniprot ID	Q8NAU1
Expiration Date	6 months from date of receipt.