

Product Datasheet

Recombinant Human IL18R1 (C-6His) (orb1906559)

Description	Recombinant Human Interleukin-18 Receptor 1 is produced by our Mammalian expression system and the t...
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-6His
Note	For research use only
Uniprot ID	Q13478
Expiration Date	6 months from date of receipt.