

Product Datasheet

Recombinant Cynomolgus PVRIG (C-Fc) (orb1906604)

Description	Recombinant Cynomolgus Transmembrane Protein PVRIG is produced by our Mammalian expression system an...
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc
Note	For research use only
Uniprot ID	A0A2K5WVV8
Expiration Date	6 months from date of receipt.