

## Product Datasheet

### Recombinant Human VPS29 Protein (orb1906746)

<b>Description</b>	Recombinant Human VPS29 Protein is produced by E. coli. The target gene encoding N68-P182 is expressed...
<b>Form/Appearance</b>	Supplied as a 0.2 µm filtered solution of PBS, 300mM NaCl, 10% Glycerol, pH 7.4
<b>Storage</b>	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Tag</b>	C-FLAG
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q9UBQ0</b>
<b>Expiration Date</b>	6 months from date of receipt.