

## Product Datasheet

# Recombinant Human Ubiquitin-conjugating enzyme E2 R2/UBE2R2 (orb1906771)

<b>Description</b>	Recombinant Human Ubiquitin-conjugating enzyme E2 R2/UBE2R2 is produced by E. coli. The target gene e...
<b>Form/Appearance</b>	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH 7.4
<b>Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Tag</b>	C-FLAG
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q712K3</b>
<b>Expiration Date</b>	6 months from date of receipt.