

Product Datasheet

Recombinant Human Protein arginine N-methyltransferase 1/PRMT1 (orb1906823)

Description	Recombinant Human Protein arginine N-methyltransferase 1/PRMT1 is produced by E. coli. The target ge...
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Uniprot ID	Q99873
Expiration Date	6 months from date of receipt.