

## Product Datasheet

# Recombinant Human Ubiquitin carboxyl-terminal hydrolase 10/USP10 (orb1906846)

<b>Description</b>	Recombinant Human Ubiquitin carboxyl-terminal hydrolase 10/USP10 is produced by E. coli. The target ...
<b>Form/Appearance</b>	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH 7.4
<b>Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Tag</b>	C-His
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q14694</a>
<b>Expiration Date</b>	6 months from date of receipt.