

Product Datasheet

Recombinant Human GAPDHS Protein (orb1906942)

Description	Recombinant Human GAPDHS Protein is produced by E. coli. The target gene encoding P20-A160 is expressed with a His tag at the C-terminus.
Form/Appearance	Supplied as a 0.2 µm filtered solution of PBS, 300mM NaCl, 10% Glycerol, pH 7.4
Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Uniprot ID	O14556
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)