

Product Datasheet

Recombinant Human SSU72 Protein (orb1906954)

| | |
|------------------------|--|
| Description | Recombinant Human SSU72 Protein is produced by E. coli. The target gene encoding P2-Y194 is expressed ... |
| Form/Appearance | Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10% Glycerol, pH 7.4 |
| Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Tag | C-His |
| Note | For research use only |
| Uniprot ID | Q9NP77 |
| Expiration Date | 6 months from date of receipt. |