

Product Datasheet

Recombinant Human Ubiquitin/ISG15-Conjugating Enzyme E2 L6/UBE2L6//UbcH8 (C-6His) (orb1906961)

Description	Recombinant Human Ubiquitin-Conjugating Enzyme E2 L6 is produced by our E.coli expression system and...
Form/Appearance	Supplied as a 0.2 um filtered solution of 50mM HEPES, 100mM NaCl, pH 8.0.
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Uniprot ID	O14933
Expiration Date	6 months from date of receipt.