

Product Datasheet

Recombinant Human Myelin P2 protein/PMP2 (orb1906990)

Description	Recombinant Human Myelin P2 protein/PMP2 is produced by E. coli. The target gene encoding S2-V132 is e...
Form/Appearance	Supplied as a 0.2 µm filtered solution of PBS, 300mM NaCl, 10% Glycerol, pH 7.4
Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Uniprot ID	P02689
Expiration Date	6 months from date of receipt.