

## Product Datasheet

# Recombinant Human Ubiquitin carboxyl-terminal hydrolase 45 (orb1907033)

|                        |  |
|------------------------|--|
| <b>Description</b>     | MADE TO ORDER   Recombinant Human Ubiquitin carboxyl-terminal hydrolase 45 is produced by E.coli exp...  |
| <b>Form/Appearance</b> | Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH 7.4   |
| <b>Storage</b>         | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles. |
| <b>Tag</b>             | C-His  |
| <b>Note</b>            | For research use only  |
| <b>Uniprot ID</b>      | <b>Q70EL2</b>  |
| <b>Expiration Date</b> | 6 months from date of receipt.   |