

Product Datasheet

Recombinant Human Methyl-CpG-Binding Protein 2/MECP2 (N-Flag&C-6His) (orb1907137)

Description	Recombinant Human Methyl-CpG-binding protein 2 is produced by our Mammalian expression system and th...
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM Tris-HCl, 300mM NaCl, 5% Trehalose, 2mM DTT, pH 8.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-Flag&C-His
Note	For research use only
Uniprot ID	P51608
Expiration Date	6 months from date of receipt.