

## Product Datasheet

# Recombinant Human Ubiquitin-Conjugating Enzyme E2 B/UBE2B/HR6B (C-6His) (orb1907205)

<b>Description</b>	Recombinant Human Ubiquitin-Conjugating Enzyme E2 B is produced by our Mammalian expression system a...
<b>Form/Appearance</b>	Supplied as a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4
<b>Storage</b>	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
<b>Tag</b>	C-His
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P63146</b>
<b>Expiration Date</b>	6 months from date of receipt.