

## Product Datasheet

# Recombinant Human Mortality Factor 4-Like Protein 2/MORF4L2/MRGX (C-6His) (orb1907208)

<b>Description</b>	Recombinant Human Mortality Factor 4-Like Protein 2 is produced by our E.coli expression system and ...
<b>Form/Appearance</b>	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
<b>Storage</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
<b>Tag</b>	C-His
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q15014</a>
<b>Expiration Date</b>	6 months from date of receipt.