

Product Datasheet

Recombinant Human HE4/WFDC2 (C-His) (orb1907216)

Description	Recombinant Human WAP Four-Disulfide Core Domain Protein 2 is produced by our Mammalian expression s...
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.2
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-6His
Note	For research use only
Uniprot ID	Q14508
Expiration Date	6 months from date of receipt.