

Product Datasheet

Recombinant Human Kallikrein 5/KLK5 (C-His) (orb1907227)

Description	Recombinant Human Kallikrein 5 is produced by our Mammalian expression system and the target gene en...
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM MES, 150mM NaCl, 10% Glycerol, pH 5.5
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-6His
Note	For research use only
Uniprot ID	Q9Y337
Expiration Date	6 months from date of receipt.