

Product Datasheet

Thalidomide-Piperazine 5-fluoride (orb1909084)

Description	Thalidomide-Piperazine 5-fluoride is a compound synthesized as an E3 ligase ligand-linker conjugate, incorporating a cereblon ligand derived from Thalidomide and utilizing a linker commonly employed in PROTAC technology.
Storage	-20°C
Note	For research use only
Formula	$C_{17}H_{17}FN_4O_4$
Purity	98.00%
MW	360.34
CAS Number	2222114-22-7
SMILES	<chem>O=C1N(C(=O)C=2C1=CC(=C(F)C2)N3CCNCC3)C4C(=O)NC(=O)CC4</chem>
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)