

Product Datasheet

Human Tumor Necrosis Factor Receptor Superfamily Member 9 (TNFRSF9) ELISA Kit (orb1909605)

Description	Human Tumor Necrosis Factor Receptor Superfamily Member 9 (TNFRSF9) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q07011
Sensitivity	32.5pg/mL
Expiration Date	Please enquire.