

## Product Datasheet

# Human Cytokine Receptor Like Factor 2 (CRLF2) ELISA Kit (orb1909622)

<b>Description</b>	Human Cytokine Receptor Like Factor 2 (CRLF2) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/mL
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma, tissue homogenates and other biological fluids.
<b>Uniprot ID</b>	<b>Q9HC73</b>
<b>Sensitivity</b>	0.052ng/mL
<b>Expiration Date</b>	Please enquire.