

Product Datasheet

Human Cyclic ADP Ribose Hydrolase (cADPRH) ELISA Kit (orb1909628)

Description	Human Cyclic ADP Ribose Hydrolase (cADPRH) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P28907
Sensitivity	0.068ng/mL
Expiration Date	Please enquire.