

Product Datasheet

Human Tumor Necrosis Factor Receptor Superfamily Member 9 (TNFRSF9) ELISA Kit (orb1909636)

Description	Human Tumor Necrosis Factor Receptor Superfamily Member 9 (TNFRSF9) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q07011
Sensitivity	32.5pg/mL
Expiration Date	Please enquire.