

Product Datasheet

Human Cytokine Receptor Like Factor 2 (CRLF2) ELISA Kit (orb1909653)

Description	Human Cytokine Receptor Like Factor 2 (CRLF2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9HC73
Sensitivity	0.052ng/mL
Expiration Date	Please enquire.