

Product Datasheet

Mouse R-Spondin 3 (RSPO3) ELISA Kit (orb1909675)

Description	Mouse R-Spondin 3 (RSPO3) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q2TJ95
Sensitivity	0.061ng/mL
Expiration Date	Please enquire.