

Product Datasheet

Human Salusin Alpha (SALa) ELISA Kit (orb1909684)

Description	Human Salusin Alpha (SALa) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	4.47pg/mL
Expiration Date	Please enquire.