
Product Datasheet

**Mouse Her2 / ErbB2 Protein, Fc Tag (MALS
verified) (orb1909725)**

Description

Mouse Her2 Protein, Fc Tag is

(HEK293). It contains AA Thr 23 - Thr 6...

Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Her2

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 96.1 kDa . The protein migrates as 120-130 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Mouse Her2 Protein, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Thr 23 - Thr 653 (Accession # P70424).

MW

96.1 kDa

Uniprot ID

P70424

Expiration Date

6 months from date of receipt.