

---

## Product Datasheet

### Cynomolgus CD24 Protein, His Tag (orb1909729)

**Description**

Cynomolgus CD24 Protein, His Tag is

(HEK293). It contains AA Asn 27 -...

**Reactivity**

Monkey

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

CD24

**Preservatives**

Powder

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

**Tag**

C-10×His

**Note**

For research use only

**Application notes**

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 5.1 kDa. The protein migrates as 27-35 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Purity**

90%

**Source**

Cynomolgus CD24 Protein, His Tag is expressed from human 293 cells (HEK293). It contains AA Asn 27 - Gly 57 (Accession # A0A7N9D2D9-1).

**MW**

5.1 kDa

**Uniprot ID****A0A7N9D2D9****Expiration Date**

6 months from date of receipt.