
Product Datasheet

**HSV-2 (strain 333) Glycoprotein B (gB), His Tag
(MALS verified) (orb1909736)**

Description

HSV-2 (strain 333) Glycoprotein B

	human 293 cells (HEK293). It conta...
Reactivity	Virus
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Glycoprotein B / gB
Preservatives	Powder
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 74.3 kDa. The protein migrates as 70-90 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	HSV-2 (strain 333) Glycoprotein B (gB), His Tag is expressed from human 293 cells (HEK293). It contains AA Glu 98 - Ala 730 (Accession # P06763).
MW	74.3 kDa
Uniprot ID	P06763
Expiration Date	6 months from date of receipt.