

---

## Product Datasheet

**Human Clec4a Protein, Fc Tag (MALS verified)  
(orb1909739)**

**Description**

Human Clec4a Protein, Fc Tag is

(HEK293). It contains AA Gln 70 - Leu...

**Reactivity**

Human

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

Clec4a

**Preservatives**

Powder

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

**Tag**

C-Fc

**Note**

For research use only

**Application notes**

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 46.4 kDa. The protein migrates as 50-55 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Purity**

95%

**Source**

Human Clec4a Protein, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Gln 70 - Leu 237 (Accession # Q9UMR7-1).

**MW**

46.4 kDa

**Uniprot ID**[Q9UMR7](#)**Expiration Date**

6 months from date of receipt.