
Product Datasheet

**Mouse Clec2e Protein, Fc Tag (MALS verified)
(orb1909741)**

Description

Mouse Clec2e Protein, Fc Tag is

(HEK293). It contains AA Val 60 - Gly...

Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Clec2e

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 43.0 kDa. The protein migrates as 53-63 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Mouse Clec2e Protein, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Val 60 - Gly 204 (Accession # Q80XD9-1).

MW

43.0 kDa

Uniprot ID

Q80XD9

Expiration Date

6 months from date of receipt.