
Product Datasheet

Cynomolgus ASGR1 Protein, Fc Tag (orb1909746)

Description

Cynomolgus ASGR1 Protein, Fc Tag is

(HEK293). It contains AA Ser 61 -...

Reactivity

Monkey

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

ASGR1

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

N-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the N-terminus. The protein has a calculated MW of 52.9 kDa. The protein migrates as 55 kDa and 60-65 when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Cynomolgus ASGR1 Protein, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Ser 61 - Leu 291 (Accession # XP_005582755.1).

MW

52.9 kDa

Uniprot ID

XP_005582755.1

Expiration Date

6 months from date of receipt.