

Product Datasheet

Interferon receptor agonist (orb1909838)

| | |
|------------------------|--|
| Description | Interferon receptor agonist is an interferon (IFN) receptor agonist. |
| Note | For research use only |
| Formula | $C_{12}H_{22}N_4O_2$ |
| MW | 254.33 |
| CAS Number | [2215120-36-6] |
| SMILES | <chem>CCCC[C@H](CO)NC1=C(OC)C=NC(NC)=N1</chem> |
| Expiration Date | 12 months from date of receipt. |