

Product Datasheet

Thalidomide-NH-PEG4-Ms (orb1909948)

Description	Thalidomide-NH-PEG4-Ms is an E3 ligase ligand-linker conjugate that incorporates Thalidomide based c...
Note	For research use only
Formula	$C_{22}H_{29}N_3O_{10}S$
MW	527.54
CAS Number	[2140807-24-3]
SMILES	<chem>O=C1N(C(CC2)C(NC2=O)=O)C(C3=C1C=CC=C3NCCOCCOCCOCCOS(C)(=O)=O)=O</chem>
Expiration Date	12 months from date of receipt.