

## Product Datasheet

### Thalidomide-NH-C2-PEG3-OH (orb1909949)

<b>Description</b>	Thalidomide-NH-C2-PEG3-OH is an E3 ligase ligand-linker conjugate that incorporates Thalidomide base...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{21}H_{27}N_3O_8$
<b>MW</b>	449.45
<b>CAS Number</b>	[2140807-23-2]
<b>SMILES</b>	<chem>O=C1N(C(CC2)C(NC2=O)=O)C(C3=C1C=CC=C3NCCOCCOCCOCCO)=O</chem>
<b>Expiration Date</b>	12 months from date of receipt.