

Product Datasheet

CP-226269 (orb1910508)

Description	CP 226269 is a potent dopamine D4 receptor agonist that induces calcium flux with EC50 of 32.0 nM. C...
Note	For research use only
Formula	$C_{18}H_{19}N_4F$
MW	310.36866
CAS Number	[220941-93-5]
SMILES	<chem>C1=CC=NC(=C1)N2CCN(CC2)CC3=CC4=CC(=CC=C4N3)F</chem>
Expiration Date	12 months from date of receipt.