

## Product Datasheet

### SIRT5 inhibitor 5 (orb1910515)

<b>Description</b>	SIRT5 inhibitor 5 is a potent SIRT5 inhibitor with an IC <sub>50</sub> value of 0.21 μM.
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>21</sub> H <sub>14</sub> ClN <sub>3</sub> O <sub>3</sub> S
<b>MW</b>	423.87
<b>SMILES</b>	<chem>CC1=NN(C/C1=C/C2=CC=C(C=C2)C(O)=O)=O)C3=NC(C4=CC=C(C=C4)Cl)=CS3</chem>
<b>Expiration Date</b>	12 months from date of receipt.